


POCO | **Getting Started**



Usage: Poco is designed to mimic two fingers so you can easily bend & adapt it to take the shape that works best for you and deliver targeted vibrations, precisely where you need it.



Charging: Connect the USB wire to the port on the back of your Poco and the other end to any USB plug for power. The red light under the power icon will blink to indicate charging. When the light turns off, charging is complete. Charge for 60 minutes before use.



App: Download the free MysteryVibe App from the Apple App Store or Google Play. Connect and use the App as a remote control for your Poco. Download different vibration patterns or create your own.



Warranty Information: Register your Poco at mysteryvibe.com/register to activate the 2 year warranty and receive new product information and offers. Discover more at mysteryvibe.com.

FDA-registered Class II therapeutic genital vibrator (21 CFR § 884.5960, Product Code KXQ). Not intended to diagnose, treat, cure, or prevent any disease or condition.

6903-2008-105-1.2